

**Clinical Laboratory**

<b>URINARY ALDOSTERONE</b>	<b>URINARY POTASSIUM</b>
<b>URINARY CHLORIDE</b>	<b>URINARY COPPER</b>
<b>URINARY CORTISOL</b>	<b>URINARY SODIUM</b>
<b>URINARY MAGNESIUM</b>	<b>URINARY UREA</b>
<b>URINARY CALCIUM</b>	<b>URINARY URIC ACID (URICURIA)</b>
<b>CREATININE CLEARANCE</b>	<b>PROTEINURIA</b>
	<b>URINARY ZINC</b>

**24-hour urine collection procedure**

Discard the first morning urine (at 7:00 AM – start of collection) and collect all subsequent urine, including the first urine of the following morning (at 7:00 AM – end of collection).

During collection, avoid exposing the container to light and temperatures above 10°C.

Only use the designated containers available from pharmacies.

For CREATININE CLEARANCE testing, a serum creatinine test (blood sample) is also required.

Milan, 13/11/2024